

Regevan[®] fill
volume **MAX**


















Regevan[®] fill
volume **PRO**

Regevan[®] fill
volume **UP**




Tentech Item

Regevan[®] fill volume UP Regevan[®] fill volume PRO Regevan[®] fill volume MAX

									
<p>10 THERMA</p> <p>Monopolar RF</p>	<p>10 THERA</p> <p>2 Line HIFU</p>	<p>10 PL</p> <p>168 Cell Fractional IPL</p>	<p>10 CELL</p> <p>3cm x 3cm CO2 Fractional Laser</p>	<p>ONE THERA</p> <p>Dual Mode Ultrasound System</p>	<p>TEN HI</p> <p>Fractional & Thermal RF</p>	<p>10 SONO</p> <p>DSS & NBS Dual Mode System</p>	<p>LIPO THERA</p> <p>Dual Mode Ultrasound System</p>	<p>10 LUX</p> <p>High Intensity & Multiple Light Therapy Device</p>	<p>10 LUX HAIR</p> <p>PhotoBioModulation Device</p>
			 <p>Rhino fill Premium</p> <p>Regevan fill volume MAX</p> <p>Regevan fill volume PRO</p> <p>Regevan fill volume UP</p> <p>Regevan fill FINE</p> <p>Monophasic Filler [for Export]</p>			<p>Regevan[®] Pn</p> <p>Polynucleotide Sodium</p>			
<p>Regevan[®] LED Mask</p> <p>Dr. Oracle LED Mask 4,474 LED Lights</p>	<p>10 mono Thread</p> <p>Mono Thread [Registered for NMPA]</p>	<p>10 cog Thread</p> <p>Asymmetrical Cog Thread</p>					<p>Regevan[®] C</p> <p>Growth Factor Ointment</p>	<p>Regevan[®] S</p> <p>Growth Factor Spray</p>	<p>Regevan[®] ALO CARE</p> <p>Growth Factor Ampoule</p>

1. High-quality Hyaluronic Acid Raw Material


Is made of High-quality Hyaluronic Acid Raw Material, listed in the U.S. Food and Drug Administration's (FDA) Drug Master Files (DMFs) for drug raw material and certified by the European Directorate for the Quality of Medicines & HealthCare (EDQM).



* Hyaluronic Acid Raw Material is produced from Non-GMO Streptococcus zooepidemicus, and animal culture medium and enzymes are not used in the entire manufacturing process.

2. Soft and Elastic Property

Patients can expect natural treatment effects from soft and highly elastic properties achieved through high G' values and tangent values. (According to Tentech's equipment measurement standards)




* When the tangent value is bigger than 1, the property is closer to viscous liquid and when it is less than 1, it is closer to elastic solid.

3. Even and Stable Injection Force




Stable injection force enables to make even and predictable treatment.

Injection Force Test Comparison



Specification

Regevan® fill volume UP Regevan® fill volume PRO Regevan® fill volume MAX

Product	Regevan® fill volume UP	Regevan® fill volume PRO	Regevan® fill volume MAX
Image			
Category	Medical Device		
HA Concentration	20mg/mL		
Viscosity	★★★	★★☆	★★★
Volume	1mL X 1Syringe	1mL X 2Syringe	
Needle	27G * 60mm * 2ea	27G * 1/2'(13mm) UTW	
Lidocaine	0.15%	0.3%	
Cross-link Completion	BDDE		
Indication	<p>By injecting cross-linked Hyaluronic acid into the dermis, it temporarily improves facial wrinkles in adults through physical recovery.</p> <p style="text-align: center;">+</p> <p>It temporarily improves mid-face volume loss in adults with moderate to severe diseases or higher.</p>		

Regevan[®] fill **Regevan[®] fill** **Regevan[®] fill** volume **UP** volume **PRO** volume **MAX**

All information in this document is subject to copyrights and other intellectual property rights of Tentech Inc. and other licensees. No modification, reproduction or copy of any or any part of this material may be made without prior written permission from Tentech Inc. and licensees. Additional information is available upon request. The following are trademarks or registered trademarks of Tentech, Inc. in Korea and other countries.

